

# Technical data ImageMaster® MF 200

	ImageMaster® MF 200
Accuracy (MTF on-axis and off-axis)	±2 % MTF on-axis ±3 % MTF off-axis
Repeatability (MTF on-axis and off-axis)	±1 % MTF on-axis ±1,5 % MTF off-axis
Optical set up	Infinite imaging set up
Sample EFL (Effective Focal Length)	18 mm ... 90 mm
Sample HFOV	6° ... 45°
Sample FOV	12° ... 90°
Dimensions (h x w x d)	
Weight	
Type	Stand alone